





Medical sealers





The utmost evolution in sterilization


The key point is to clearly define the process critical parameters. In order to guarantee a seal which is conform to EN 868-5 standards, the sealer must be able to monitor and record the temperature, the force applied and the time

cycle of the seal process. Such an equipment can be part of a sterile packaging process which is validatable according to EN ISO 11607-2:2006 standards. The wide Gandus Saldatrici range, from the brand new Pura to the extraordinary for its unicity VSeal.it, from the rotary sealers minirò serie, to the impulse sealers Medical H and SI series, all guarantee a seal which fully complies with the EN ISO 11607-2:2006 standard.

Our vast experience, and ongoing commitment to ensure that our wide range of equipment meets the latest requirements, allows us to offer not only a range of products which can assure the highest level of safety and service, but also a range of services totally dedicated to the Validation Process.

Safety is our DNA.

Our identity, a tradition based on quality design and tailor made manufacturing, all Italian



Few people know that Gandus history goes back to 1885, when a distinguished gentleman named Prospero Gandus founded his company in Milan. The production began with paper bags to contain the needles for the sewing machines. Craftsmanship and precision, all Italian. Since its foundation Gandus has been a manufacturing company. And since its foundation Gandus has been managed by Prospero gandus first and his descendance then.

In 1954, Gandus began to enter the world of packaging machinery production, adding the word "Saldatrici" (in English Sealers) to the brand and affirming itself immediately on the market thanks to professionalism, reliability, high quality constantly in time.

Gandus Saldatrici has today 3 specialized divisions: for Bags and tubes sealers production, for Sterilization pouches sealers production, for Automatic packaging systems production

All our machines are designed in-house, making use of the latest 3D design software, all the control electronics of our machines and software operation are designed and developed directly from Gandus, the know-how of the operation of our machines is resident in Gandus. Our mission is to design and manufacture top quality machines, to achieve this even the smallest mechanical component must be made with the utmost precision. Our technicians assemble the sealers that you use every day, working with the care of a master watchmaker, assembling hundreds of parts and components to get wonderful concentrates of high technology.

Our identity is being a team, take our responsibilities, have respect for everyone, communicate aiming to getting better and better, work professionally.

Our success is based on entrepreneurship, entrepreneurial skills, commitment, continuity and professional fairness.

This is our world, this is our way of conceiving design and production of the equipments we install all over the world.





Packaging Process Validation

The whole sterilization process is today fully regulated, to get to the highest levels in quality and safety. Packaging is a key element in the Sterilization process. Sealing has to ensure a reliable barrier against microorganism. A non conforming seal invalidates the whole Sterilization process.

Only sealing by use of Validatable ISO 11607-2 sealing devices, guarantees a perfect seal tightness. Only machines able to control and monitor critical parameters can be part of a Validatable sterile packaging process according to standard EN ISO 11607-2.

Gandus Saldatrici range of products guarantees the strict observance of a qualitative standard to assure a consistant sealing integrity and the repeatability of the process, Validatable EN ISO 11607-2.

Safety

Our mission is taking care of the whole packaging process to guarantee the highest levels of quality able to comply to the rules requirements in terms of safety in sterilization.

From the kind of seal to the parameters controls, from the routine daily tests to verify the system of closure according to EN ISO 11607,1.2 and ASTM F1929, to the accessories related to Traceability and to the services dedicated to Validation Process.

Sealers have to assure the final result is steady and error free.

Communication Interfaces

Gandus sealers are equipped with dedicated interface hardwares to allow links with external communication systems to create Validation Process editing.

Traceability

Traceability is the ability to identify, monitor and recall the history of all medical devices that come into contact with the patient or that fall into the diagnostic and therapeutic process involved. The traceability system has the primary purpose of protecting the patient's safety and can be manual, semi-computerized or fully computerized.

SecureScan Data Upload Software

Data upload software to allow you to create, print and download bar codes with all the sealing parameters of each ID pouch on an external PC. Simple and extremely safe data upload software able to read the customized packaging data and automatically set the suitable machine function and parameters.

Securscan Data Upload Software allows you fully operate it with an extremely safe and simple data upload software, able to read the customized packaging data and automatically set the suitable machine function and parameters. The system is extremely easy to install and can work on any PC. It allows a very quick and 100% safe and correct setting of your sealing device.





Packaging

In relation to the sterilization methods and to the kind of medical device, it is necessary to choose the most suitable packaging to guarantee a perfect sterility of the content.

Gandus Saldatrici offers products able to seal all kind of material, while assuring the highest levels of quality to preserve sterility.

HI Human Interface

The Packaging process is a continuous interaction between man and sealer. It is fundamental that this interaction takes place within the highest level of efficacy, and functionality, in order to ensure the utmost level of safety. Our machines are designed and manufactured offering different kind of interfaces able to meet every different requirement.

Room Saving

Our ongoing target is to design sealers that combine minimum sizes with all highest performances, integrating all functions to allow the optimum use even in reduced spaces.

Cleaning and Disinfection

All the materials used, stainless steel, anodized aluminium, together with a smooth design, allow a quick, easy and effective cleaning and disinfection.

Print Facility

Almost all Gandus medical sealers allow you to print directly on the pouch all important information related to sterile tools packaging process, such as production and expiry date, batch, operator, sealing temperature etc. All these data can be customizely printed, thanks to an advanced technology that enable the system to recognize pouches dimensions and automatically adjust character type and dimensions according to the available print area.

Energy Saving

Our machines are designed and manufactured aiming to the optimization of processes, whilst respecting the need of energy saving.

Vacuum Packaging

Hygienic, space saving and better shelf life requirements are perfectly satisfied by Gandus sealers vacuum systems. Thanks to vacuum packaging, the package size is reduced by 90% of its volume. Highly suitable to be used for storage of cloths, uniforms, bed linen, that need minimum space.



Index

· Rotary Sealers Pura	7
 Pouch Maker Vseal.it 	8
 Pouch Packaging Station ONE4ALL 	9
 Rotary Sealers minirò H-lan - minirò H-lan touch minirò H-seal Evo, minirò H-seal VP Evo minirò H-data Evo, minirò minirò H-net Evo 2.0 minirò poly VP Evo, rebi Evo, rebi Prt Evo, rebi VP Evo 	10 13 14 15 16
• Impulse Sealers Medical H 460/610 - Medical H 460/610 Digital AccuTemp	11
 Vacuum Impulse Sealers SI Med - SI Vacuummed 520/720/1020 	12
• Dental Sealers DN 350, GD 301 <i>Evo</i>	17
 Gandus Routine Tests BlueTest.it SealCheckIT 	18
· Accessories	19 - 20
 Gandus Services to Validation Process Calibration Service Peel Strength Test Validation Consultancy Service (IQ, QQ) 	21
• Technical features	22 - 26







Rotary Sealers

Pura

Product design and engineering developed with the scientific support of Politecnico di Milano, Design Department.



Accurate temperature control thanks to new generation control hardware and software. Fully complying with ISO 11607-2 standard and international guidelines EN-ISO TS 16775.

Safety

The reliable and well known Gandus 15 mm wide multiline seal guarantees a perfect seal strength with an easy and safe peelability.

Communication Interfaces

USB and RS 232 interfaces allows communication with several external devices.

Traceability

Pura allows fully traceability of all sealing parameters on a standard USB stick.

SecureScan Data Upload Software

Operating with a scanner, Securscan Data Upload Software allows to upload all data on the sealer.

Packaging

Pura is suitable for the sealing of all known, sealable, pre-made sterile barrier systems, uncoated Tyvek®, as well as the newest generation of packaging material.

HI - Human Interface

4,3" Wide Color touch screen and easy to use software are the perfect Human Interface.

Room Saving

Just 0,15 square meter are required to place a Pura on a packaging table.

Cleaning and Disinfection

Stainless steel, anodized aluminium, together with an unbelievable ergonomic design, allows and easy and effective cleaning and disinfection.

Print Facility

Last generation 9 dot matrix printer enables to print on packaging of all the required data.

Automatic adjustment of font size to match pouch width.

Energy saving















Vseal.it

Pouch maker



Packaging Process Validation

Accurate temperature control thanks to new generation control hardware and software. Fully complying with ISO 11607-2 standard and international guidelines EN-ISO TS 16775.

Safety

The Vseal.it is the unique pouch maker that allows manifacture sterilisable medical pouches with Chevron Seal from standard reels. Chevron seal means peelability, safety, correct pouch opening and aseptic presentation.

Communication Interfaces

USB and RS 232 interfaces allows communication with several external devices.

Traceability

Vseal.it allows full traceability of all sealing parameters on a standard USB stick.

SecureScan Data Upload Software

Operating with a scanner, Securscan Data Upload Software allows to upload all data on the sealer.

Packaging

Vseal.it is suitable for manufacturing BOPP or Tyvek® pouches, as well as the newest generation of packaging material.

HI - Human Interface

4,3" Wide Color touch screen and easy to use software are the perfect Human Interface.

Room Saving

Using just one Vseal.it it is possible to save hundreds square meters of warehousing space to stock your sterilization packaging pouches.

Cleaning and Disinfection

Stainless steel top cover allows an easy and effective cleaning and disinfection.

Energy saving













Packaging Station

ONE4ALL



Packaging Process Validation

ONE4ALL is a complete Pouch Packaging Station designed to gather together in one unique element all necessary devices for an easy and perfect Validated Packaging Process.

Safety

Its smooth design and its utmost ergonomics offer the utmost safety for operators.

Communication Interfaces

ONE4ALL offers both electrical and Ethernet connections.

Traceability

The Ethernet interface allows communication with dedicated softwares for a total traceability of all required data.

Packaging

Designed for a perfect use and storage of all required packagings: pouches, big pouches, reels. Ideal for BOPP, Tyvek®, SMS, plastics, header bags.

HI - Human Interface

Buttons to electrically adjust the height according to operator's height. A unique electric panel allows an easy connection of all equipments.

Room Saving

It deals with a compact block 2mts width x 2 mts length designed for a perfect use of all 4 devices, offering wide workspace. No area limits, it can be placed anywhere as it moves on wheels, specific for hospital use.

Ergonomics

Designed thinking about the need of an easy lodging of all required materials and packed products.

Electrically adjustable in height according to operator's height. No electric connection limits as all equipments are connected to a unique electric panel.

Cleaning and Disinfection

The use of alluminium and stainless steel allows an easy and effective cleaning and disinfection.

Energy saving











Rotary Sealers

minirò H-lan touch



Packaging Process Validation

Microprocessor temperature controller to ensure an accurate temperature control. Fully complying with ISO 11607-2 standard and international guidelines EN-ISO TS 16775.

Safety

The reliable and well know Gandus 15 mm wide multiline seal guarantee a perfect seal strength with an easily and safety peelability.

Communication Interfaces

Ethernet, USB and RS 232 multiple interfaces allow communication with several external devices and remote control of the sealer.

Traceability

minirò H-lan allows fully traceability of all sealing parameters on board or an external computer server.

Packaging

The minirò H-lan is suitable for the sealing of all known, sealable, pre-made sterile barrier systems, uncoated Tyvek®, as well as the newest generation of packaging material.

HI - Human Interface

7" Wide Color touch screen and easy to use software are the widest Human Interface available on the market of medical rotary sealers.

Room Saving

With only $62\,\mathrm{cm}$ width, the minirò H-lan is equipped with a built-in double line printer.

Cleaning and Disinfection

Stainless steel, together with a unique design, allow an easy and effective cleaning and disinfection.

Print Facility

18 dots matrix printer to print 1 line, 2 lines or barcode on packaging.

Energy saving









Impulse Sealers

Medical H 460/610 Medical H 460/610 DIGITAL ACCUTEMP

Packaging Process Validation

Highest world technology used for temperature controller of impulse sealers ensures the most accurate temperature control. Fully complying with ISO 11607-2 standard and international guidelines EN-ISO TS 16775.

Safety

8 mm wide flat seal guarantees a perfect seal strength with an easy, safe and effective peelability.

Communication Interfaces

RS 232 interface allows communication with a dedicated software for traceability.

Traceability

A dedicate traceability software is available for Gandus Medical Digital Accutemp.

Packaging

The Medical Digital Accutemp is suitable for the sealing of all known, sealable, pre-made sterile barrier systems, uncoated Tyvek®, as well as the newest generation of packaging material.


HI - Human Interface

Digital display 20 characters, 4 lines and easy to use software allows a user friendly approach to our sealer.

Cleaning and Disinfection

Stainless steel top cover allows an easy and effective cleaning and disinfection.

Energy saving
















Vacuum Impulse Sealers

SI Med SI Vacuummed 520/720/1020



Highest world technology from temperature controller of impulse sealers ensures the most accurate temperature control. Fully complying with ISO 11607-2 standard and international guidelines EN-ISO TS 16775.

Safety

8 mm wide flat seal guarantee a perfect seal strength with an easily and safety peelability.

Communication Interfaces

A DIAG interface allows communication with a dedicated software for traceability.

Traceability

A dedicate traceability software is available for our GANDUS SI serie.

Vacuum Packaging

SI Vacuum sealer employs a powerful vacuum nozzle with a double sealing bars to ensure that products are packed in airless or modified atmosphere environment. SI Vacuum sealer is able to guarantee up to 70% vacuum level.

HI - Human Interface

Digital display 20 characters, 4 lines and easy to use software allows a user friendly approach to our sealer.

Cleaning and Disinfection

Fully stainless steel allows an easy and effective cleaning and disinfection.

Energy saving














Rotary Sealers

minirò H-seal Evo minirò H-seal VP Evo



The model minirò H-seal VP *Evo*, fully complies with ISO 11607-2 standard and international guidelines EN-ISO TS 16775.

Safety

The reliable and well known Gandus 15 mm wide multiline seal guarantees a perfect seal strength with an easy and safe peelability.

Communication Interfaces

A USB interface allows connection to USB pendrive standard.

Traceability

The minirò H-seal VP *Evo* allows fully traceability of all sealing parameters on a standard USB stick.

Packaging

All miniro H-seal *Evo* models are suitable for the sealing of all known, sealable, pre-made sterile barrier systems, uncoated Tyvek®, as well as the newest generation of packaging material.

HI - Human Interface

User friendly 16 characters on 2 lines digital display and a wide and clear keyboard.


Cleaning and Disinfection

Stainless steel top cover allows an easy and effective cleaning and disinfection.

Print Facility

On the model minirò H-seal VP *Evo* all the requested data could be printed with a label printer linked to the sealer.

Energy saving











Rotary Sealers

minirò H-net Evo 2.0 minirò H-data Evo

Packaging Process Validation

The model minirò H-net *Evo* 2.0, Fully complies with ISO 11607-2 standard and international guidelines EN-ISO TS 16775.

Safety

The reliable and well known Gandus 15 mm wide multiline seal guarantees a perfect seal strength with easy and safe peelability.

Communication Interfaces

On the model minirò H-net *Evo* 2.0 USB and RS 232 multiple interfaces allow communication with several external devices.

Traceability

The model minirò H-net *Evo* 2.0 allows fully traceability of all sealing parameters on a standard USB stick.

SecureScan Data Upload Software

Operating with a scanner, Securscan Data Upload Software allows to upload all data on the sealer.

Packaging

These machines are suitable for the sealing of all known, sealable, pre-made sterile barrier systems, uncoated Tyvek®, as well as the newest generation of packaging material.

HI - Human Interface

The model minirò H-net *Evo* 2.0 is equipped with an user friendly 20 characters on 4 lines digital display and a wide and clear keyboard. The model minirò H-data *Evo* is equipped with the well known 16 characters on 2 lines digital display and a wide and clear keyboard.

Cleaning and Disinfection

Stainless steel top cover allows an easy and effective cleaning and disinfection.

Print Facility

Last generation 9 dot matrix printer allows to print on packaging all the requested data.

Energy saving













Rotary Sealer

minirò H-poly VP Evo



Packaging Process Validation

The model minirò H-poly VP *Evo*, fully complies with ISO 11607-2 standard and international guidelines EN-ISO TS 16775.

Safety

A 9 mm wide knurled seal guarantees a perfect seal strength with easy and safe peelability.

Communication Interfaces

A USB interface allows connection to USB pendrive standard.

Traceability

The model minirò H-poly VP *Evo* allows fully traceability of all sealing parameters on a standard USB stick.

Packaging

The model minirò H-poly VP *Evo* is suitable for the sealing LDPE pouches, header bags and other all known, sealable, pre-made sterile barrier systems, uncoated Tyvek®.

HI - Human Interface

The model minirò H-poly VP *Evo* is equipped with the well-known 16 characters on 2 lines digital display and a wide and clear keyboard.

Cleaning and Disinfection

Stainless steel top cover allows an easy and effective cleaning and disinfection.

Print Facility

All the requested data can be printed with a label printer linked to the sealer.

Energy saving









Rotary sealers

rebi Evo - rebi Prt Evo rebi VP Evo

Packaging Process Validation

Gandus model rebi VP *Evo*, fully complies with ISO 11607-2 standard and international guidelines EN-ISO TS 16775.

Safety

The reliable and well known Gandus 12 mm wide multiline seal guarantees a perfect seal strength with an easy and safe peelability.

Communication Interfaces

On Gandus rebi VP *Evo*, a USB interface allows connection to USB pendrive standard.

Traceability

Gandus rebi VP *Evo* allows fully traceability of all sealing parameters on a standard USB stick.

Packaging

All Gandus rebi models are suitable for the sealing of all known, sealable, pre-made sterile barrier systems, uncoated Tyvek $^{
m ll}$, as well as the newest generation of packaging material.

HI - Human Interface

User friendly 8 characters on 2 lines digital display and a wide and clear keyboard.

Cleaning and Disinfection

Stainless steel top cover allows an easy and effective cleaning and disinfection.

Print Facility

On the model rebi Prt *Evo* a 9 dot matrix printer allows to print on packaging all the requested data.

Upon request rebi VP *Evo* model can be equipped with integrated printer facility.

Energy saving







Dental Sealers

DN 350 GD 301 Evo

Safety

 $12~\mathrm{mm}$ wide flat seal on GD $301~\mathrm{Evo}$ and $8~\mathrm{mm}$ wide flat seal on DN 350, guarantees a perfect seal strength with easy and safe peelability.

Packaging

Our dental sealers are suitable for the sealing pre-made sterile barrier systems for steam sterilizers.

Cleaning and Disinfection

The frame of the GD 301 *Evo* is a plastic molding to allow an easy and effective cleaning and disinfection.

Energy saving



Gandus Routine Tests

Gandus offers routine tests to enable you an easy, quick and safe daily control of your sealing device in order to comply **EN ISO 11607.1.2 and ASTM F1929**.

Blue Test.it

New generation indicator test for routing testing of your sealing devices. According to EN ISO 11607.1.2 and ASTM F1929:

Specially developed to check the integrity of seals produced with your sealers. Suitable for Paper and Tyvek® backed laminated pouches. It allows a safe, easy and guick ink test.

Safe - no damage of the seal because it does not use any pipette or other objects to fill the pouch, and eliminates any risk of leakage or migration of the methylene blue dye.

Easy - easy to prepare and easy to check the result.

Quick - after no longer than 5 seconds per seal you can see if the seal made with your device is perfect, secure and safe.









Seal CheckIT

Test system for the daily testing of sealing device and for Operational Qualification according to EN ISO 11607-2:2006 process Validation.

Guidelines of the WFHSS are advising to test each medical heat sealer at the beginning of the day after the heat sealer has heated up. Based upon the initial test the heat sealer can be released for the day. The Gandus Saldatrici Seal CheckIT is one way to comply to the ISO 11607 - 2 requirements for Operational Qualification (OQ) and daily testing of your Medical Heat sealers.













ACCESSORIES

	Accessories	Sealer model
	Front tray	Pura
	Front tray	rebi serie
TIL	Front tray	Medical H 460/610, Medical H ATDP 460/610
	Front tray	DN 350
100	Indipendent motorized conveyor band	minirò H-poly VP <i>Evo</i>
	PS 654 Front tray	minirò H-seal <i>Evo</i> , minirò H-seal VP <i>Evo</i> , minirò H-net, minirò H-data <i>Evo</i>
	PS 70 Front tray	minirò H-lan, minirò H-lan Touch
	Slip roll	Pura
	RL 79 Slip roll	minirò H-lan, minirò H-lan Touch
	RL 80 Slip roll	minirò H-seal <i>Evo</i> , minirò H-seal VP <i>Evo</i> , minirò H-net, minirò H-data <i>Evo</i> , rebi serie






ACCESSORIES

Accessories

Sealer model



PR H40 Roll holder 1 level, complete with cutting device

All models



PR H50 Roll holder 2 levels, complete with cutting device

All models



PRN 2000 Label printer

minirò H-seal VP *Evo*, minirò H-net, minirò H-lan, minirò H-lan Touch, minirò H-poly VP *Evo*, Pura



Securscan DataUpload Software
with Barcode scanner

minirò H-net, Pura, Vseal.it



SIS

Stand adjustable in height, to allows vertical sealing

SI Med - SI Vacuummed 520/720/1020



SIPO

Adjustable work table (only with SIS)

SI Med - SI Vacuummed 520/720/1020



SIRS

Bag stretcher (only on 1020 version)

SI Med - SI Vacuummed 1020

Gandus Services to Validation Process

To obtain and maintain safety levels as from rules requirements it is necessary that sealing process is steadily monitored, in order to prevent any deviation from set qualitative standards.

Our experience and ongoing commitment to ensure that all Gandus Saldatrici equipment meets the latest requirements, allows us to offer not only a range of products which can assure the highest level of safety as from ISO 11607-2, but also a service totally dedicated to the Validation Process on all Gandus validatable medical sealers.




Calibration Service

Periodic Inspection on Medical Sealers and subsequent Calibration Certificate issue.

As part of our after sales service we are able to provide Periodic Inspection on Medical Sealers including the issue of a Machine Calibration Certificate, checking the machine sealing parameters status and ensuring compliance with ISO 11607.2, EN 868-5 requirements and to Validation process Qualification Performance (QP). All service and/or repair work is guaranteed and only carried out by our trained engineers, by means of calibrated machinery.

Gandus can also offer a reminder system on contracted work.



Peel Strength Test

Heat sealed joint sealing tests as required by UNI EN 868-5 and Qualification Performance (QP) procedures.

The UNI EN 868-5 Directive specifies the checks to be carried out on the envelope and defines the test methods to be used on the joint heat-sealed envelopes and medical rolls. By using calibrated instrumentation, Gandus Saldatrici offers Peel Strength Test service on sealing, both on sterilized envelopes and on non-sterilized envelopes, in accordance with the requirements of UNI EN 868-5 and the Validation Process Performance Qualification (QP).

Validation Consultancy Service (IQ, QQ)

Installation and Operation Qualifications are parti integranti of the Validation Process and set qualitative benchmarks that cannot be bypassed. IQ - Process of obtaining and documenting evidence that equipment has been provided and installed in accordance with its specification. Gandus Saldatrici support customers in assuring the perfect working conditions are satisfied with the ultimate aim of the utmost safety.

OQ - Process of obtaining and documenting evidence that installed equipment operates within predetermined limits when used in accordance with its operational procedures. Gandus Saldatrici support customers in assuring the Validation required pouch tightness is reached within the whole range of sealing parameters data tolerance.

If outputs are positive, IQ and OQ Certificates are issued, together with Tightness Reports (OQ).





MODELS	V seal.it	Pura minirò H-lan touch	minirò H-lan minirò H-net Evo minirò H-seal VP Evo minirò-poly VP rebi VP Evo	minirò H-data Evo minirò H-seal Evo rebi Prt Evo rebi Evo	Medical H- 460/610 Medical H 460/610 Digital Accutemp SI Med / SI Vacuummed 520/720/1020	DN 350 GD 301 <i>Evo</i>
GENERAL FEATURES		Т	T			
Aluminium		Pura				
Stainless Steel frame	•	•	•	•	•	
Painted Steel frame	•	lan touch	•	•	Medical H	DN 350
Injection plastic moulding						GD 301 <i>Evo</i>
Built-in adjustable guide		•	•	•		
Adjustable free edge		•	•	•	•	•
Built-in cutting device	•				•	•
SEALING						
Rotary sealing		•	•	•		
Impulse sealing					•	DN 350
Constant heat	•					GD 301 <i>Evo</i>
Kind of sealing	Chevron multiline	Multiline (flat option)	Multiline (flat option)	Multiline (flat option)	Flat	Flat
Sealing width (mm)		·	poly VP knurled (flat option)			
	11,5	15	lan, net, seal VP = 15 rebi VP = 12 - poly VP = 9	data, seal = 15 rebi Prt, rebi = 12	8	DN 350 = 8 GD 301 <i>Evo</i> = 12
Sealing length (mm)	300	Unlimited	Unlimited	Unlimited	Medical H = 460/610 Medical H D.A. = 460/610 SI Med/Vmed = 520/720/1020	DN 350 = 350 GD 301 <i>Evo</i> = 300
Distance sealing-medical product (mm)		≥30 Din rules	≥30 Din rules	≥30 Din rules	≥30 Din rules	≥30 Din rules
HEAT SEALABLE MATERIALS						
Pouches/reels EN 868-5	•	•	•	•	•	•
Tyvek® pouches/reels	•	•	•	•	Medical H D.A., SI Med/Vmed	
Paper pouches EN 868-4		•	•	•	IVIEUICAI II D.A., SI IVIEU/VIIIEU	•
Gusseted pouches/reels	•	•	•	•	•	•
Laminated AL pouches		•	•	•	•	
		•	•	•	•	
PP pouches		•				DN 350
Header bags			poly VP		•	DN 350
Polyethylene			poly VP		•	DN 350
Ultra® pouches/reels	•	•	•	•	Medical H D.A., SI Med/Vmed	
ISO 11607-2						
Sealing validation through automatic	•	•	•		Medical H D.A., SI Med/Vmed	
critical sealing parameters control					ivieulcai ii D.A., Si ivieu/villeu	
Temperature control						
	•	•		•	Madical II D.A. Cl Mad//mad	
	•	•	•	•	Medical H D.A., SI Med/Vmed	
Force control	•	•	•	•	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	
Speed control				•		
Speed control Time control	•	•	•	•	Medical H D.A., SI Med/Vmed	•
Speed control Time control Sealing parameters recording	•	•	•	•	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	•
Speed control Time control Sealing parameters recording Sealing parameters visualization	•	•	•		Medical H D.A., SI Med/Vmed	•
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe	•	•	•	• seal, rebi	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	•
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check	•	•	•	seal, rebi	Medical H D.A., SI Med/Vmed	•
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check Alarm and stop temperature	•	•	•		Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	•
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check Alarm and stop temperature tolerance +/-5 °C (DIN 58953-7)	•	•	•	seal, rebi • •	Medical H D.A., SI Med/Vmed	•
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check Alarm and stop temperature tolerance +/5 °C (DIN 58953-7) Alarm and stop adjustable tolerance	•	•	•	seal, rebi	Medical H D.A., SI Med/Vmed	•
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check Alarm and stop temperature tolerance +/5 °C (DIN 58953-7) Alarm and stop adjustable tolerance Force alarm and stop	•	•	•	seal, rebi • •	Medical H D.A., SI Med/Vmed	•
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check Alarm and stop temperature tolerance +/5 °C (DIN 58953-7) Alarm and stop adjustable tolerance	•	•	•	seal, rebi • •	Medical H D.A., SI Med/Vmed	•
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check Alarm and stop temperature tolerance +/-5 °C (DIN 58953-7) Alarm and stop adjustable tolerance Force alarm and stop Speed/Time alarm and stop CERTIFICATIONS/SAFETY	•	•	•	seal, rebi • •	Medical H D.A., SI Med/Vmed	•
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check Alarm and stop temperature tolerance +/5 °C (DIN 58953-7) Alarm and stop adjustable tolerance Force alarm and stop Speed/Time alarm and stop CERTIFICATIONS/SAFETY UNI EN ISO 11607-2:2006	•	•	•	seal, rebi • seal, rebi	Medical H D.A., SI Med/Vmed	
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check Alarm and stop temperature tolerance +/5 °C (DIN 58953-7) Alarm and stop adjustable tolerance Force alarm and stop Speed/Time alarm and stop CERTIFICATIONS/SAFETY UNI EN ISO 11607-2:2006 UNI EN 868-5:2009	•	•	•	seal, rebi • •	Medical H D.A., SI Med/Vmed	•
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check Alarm and stop temperature tolerance +/5 °C (DIN 58953-7) Alarm and stop adjustable tolerance Force alarm and stop Speed/Time alarm and stop CERTIFICATIONS/SAFETY UNI EN ISO 11607-2:2006 UNI EN 868-5:2009 EN ISO TS 16775	•	•	•	seal, rebi • seal, rebi	Medical H D.A., SI Med/Vmed	•
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check Alarm and stop temperature tolerance +/5 °C (DIN 58953-7) Alarm and stop adjustable tolerance Force alarm and stop Speed/Time alarm and stop CERTIFICATIONS/SAFETY UNI EN ISO 11607-2:2006 UNI EN 868-5:2009 EN ISO TS 16775 DIN 58953-7	•	•	•	seal, rebi • seal, rebi	Medical H D.A., SI Med/Vmed	
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check Alarm and stop temperature tolerance +/-5 °C (DIN 58953-7) Alarm and stop adjustable tolerance Force alarm and stop Speed/Time alarm and stop CERTIFICATIONS/SAFETY UNI EN ISO 11607-2:2006 UNI EN 868-5:2009 EN ISO TS 16775 DIN 58953-7 Directive 2006/42/EC	•	•	•	seal, rebi • seal, rebi	Medical H D.A., SI Med/Vmed	•
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check Alarm and stop temperature tolerance +/-5 °C (DIN 58953-7) Alarm and stop adjustable tolerance Force alarm and stop Speed/Time alarm and stop CERTIFICATIONS/SAFETY UNI EN ISO 11607-2:2006 UNI EN 868-5:2009 EN ISO TS 16775 DIN 58953-7 Directive 2006/42/EC Directive 2014/30/EU	•		•	seal, rebi • seal, rebi	Medical H D.A., SI Med/Vmed	•
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check Alarm and stop temperature tolerance +/-5 °C (DIN 58953-7) Alarm and stop adjustable tolerance Force alarm and stop Speed/Time alarm and stop CERTIFICATIONS/SAFETY UNI EN ISO 11607-2:2006 UNI EN 868-5:2009 EN ISO TS 16775 DIN 58953-7 Directive 2006/42/EC Directive 2014/30/EU Directive 2014/30/EU	•		• • • • • • • • • • • • • • • • • • •	seal, rebi • seal, rebi • seal, rebi •	Medical H D.A., SI Med/Vmed	• • DN 350
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check Alarm and stop temperature tolerance +/-5 °C (DIN 58953-7) Alarm and stop adjustable tolerance Force alarm and stop Speed/Time alarm and stop CERTIFICATIONS/SAFETY UNI EN ISO 11607-2:2006 UNI EN 868-5:2009 EN ISO TS 16775 DIN 58953-7 Directive 2006/42/EC Directive 2014/30/EU Directive 2014/35/EU UNI/TR 11408:2011	•		• • • • • • • • • • • • • • • • • • •	seal, rebi • seal, rebi • seal, rebi • • • • • • •	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	• • DN 350
Speed control Time control Sealing parameters recording Sealing parameters visualization Overheating temperature probe Sealer functioning diagnostic and check Alarm and stop temperature tolerance +/-5 °C (DIN 58953-7) Alarm and stop adjustable tolerance Force alarm and stop Speed/Time alarm and stop CERTIFICATIONS/SAFETY UNI EN ISO 11607-2:2006 UNI EN 868-5:2009 EN ISO TS 16775 DIN 58953-7 Directive 2006/42/EC Directive 2014/30/EU Directive 2014/35/EU	•		• • • • • • • • • • • • • • • • • • •	seal, rebi • seal, rebi • seal, rebi • • • • • • •	Medical H D.A., SI Med/Vmed	• • DN 350



MODELS PARAMETERS	Vseal.it	Pura minirò H-lan touch	minirò H-lan minirò H-net Evo minirò H-seal VP Evo minirò-poly VP rebi VP Evo	minirò H-data Evo minirò H-seal Evo rebi Prt Evo rebi Evo	Medical H- 460/610 Medical H 460/610 Digital Accutemp SI Med / SI Vacuummed 520/720/1020	DN 350 GD 301 <i>Evo</i>
Sealing speed (mt/min)	0,4 - 7	Pura 10 lan touch 3 -13	lan, net, seal VP = 10 poly VP, rebi VP = 6	data, seal = 10 rebi = 6		
Sealing time (sec.) Max sealing temperature (°C)	230	200	200 poly VP = 220	200	Medical H D.A.=250 SI Med/Vmed=200	DN 350
Reference Sealing force (N)	75	100	lan, net, seal VP = 100 poly VP, rebi VP = 85	Self-adjustable	Medical H / H D.A. about 850 N SI Med/Vmed about 750 N	GD 301 <i>Evo</i> Self-adjustable
Temperature tolerance	± 1%	± 1%	± 1%	± 1%	Medical H D.A., SI Med/Vmed ± 1%	
Force Range (N)	55-100	80-120	lan, net, seal VP = 80 - 120 poly VP, rebi VP = 70 - 102	2170	Medical H D.A.,= 820-920 SI Med/Vmed = 700-1050 (depending on the model)	
Overheating control	•	•	•	•		
CALIBRATION Heating autocalibration function Automatic self check of thermocouple battery and	•	•	•	seal / rebi ●	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	
remote server condition at start-up Temperature calibration Force calibration	•	•	•	•	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	
Speed calibration	•	•	•		Wiedical II D.A., of Wed, vinea	
Time calibration			-		Medical H D.A., SI Med/Vmed	
CULOTIONO						
FUNCTIONS Microprocessor Multilanguage menu	•	•	•	•	Madical H D A SI Mad//mad	
Microprocessor Multilanguage menu					Medical H D.A., SI Med/Vmed	
Microprocessor Multilanguage menu Starting photocell		•	• (no poly VP)	•	-	
Microprocessor Multilanguage menu		•	•	•	-	
Microprocessor Multilanguage menu Starting photocell Printing photocell	•	•	• (no poly VP)	•	Medical H D.A., SI Med/Vmed	
Microprocessor Multilanguage menu Starting photocell Printing photocell Speed adjustment	•	• • • • Ian touch	• (no poly VP)	● • data	-	
Microprocessor Multilanguage menu Starting photocell Printing photocell Speed adjustment Temperature setting	•	lan touch	• (no poly VP)	● • data	Medical H D.A., SI Med/Vmed	DN 350
Microprocessor Multilanguage menu Starting photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting	•	lan touch	• (no poly VP)	● • data	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	DN 350
Microprocessor Multilanguage menu Starting photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting	•	lan touch	• (no poly VP) Ian, net	data	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	DN 350
Microprocessor Multilanguage menu Starting photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting	•	lan touch • • • • • • • • • • • • • •	• (no poly VP) Ian, net	data data	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	DN 350
Microprocessor Multilanguage menu Starting photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting Packaging programs setting	•	lan touch • • • • • • • • • • • • • •	• (no poly VP) lan, net • (and poly VP) lan, net	data data data data	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	DN 350
Microprocessor Multilanguage menu Starting photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting Packaging programs setting Expiry dates setting	•	lan touch o o o o o o o o o	• (no poly VP) lan, net •	data data data data	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	DN 350
Microprocessor Multilanguage menu Starting photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting Packaging programs setting Expiry dates setting Stand by	•	lan touch o o o o o o o o o o o o o	• (no poly VP) lan, net • lan, net, seal VP, poly VP lan, net, seal VP, poly VP	data data data data	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	
Microprocessor Multilanguage menu Starting photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting Packaging programs setting Expiry dates setting Stand by Energy saving	•	lan touch o o o o o o o o o o o o o	Ian, net, seal VP, poly VP	data data data data data	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	DN 350
Microprocessor Multilanguage menu Starting photocell Printing photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting Packaging programs setting Expiry dates setting Stand by Energy saving Safety thermostat	•	lan touch lan touch	Ian, net, seal VP, poly VP Ian, net, seal VP, poly VP	data data data data data	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	
Microprocessor Multilanguage menu Starting photocell Printing photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting Packaging programs setting Expiry dates setting Stand by Energy saving Safety thermostat On board battery ensures data retention even when the machine is switched off	•	lan touch lan touch	• (no poly VP) lan, net • lan, net, seal VP, poly VP lan, net, seal VP, poly VP • • • • • • • • • • • • •	data data data data data	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	
Microprocessor Multilanguage menu Starting photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting Packaging programs setting Expiry dates setting Stand by Energy saving Safety thermostat On board battery ensures data retention even when the machine is switched off Peripheries condition and diagnostic check up at start-up	•	lan touch lan touch	Ian, net, seal VP, poly VP	data data data data data	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	
Microprocessor Multilanguage menu Starting photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting Packaging programs setting Expiry dates setting Stand by Energy saving Safety thermostat On board battery ensures data retention even when the machine is switched off Peripheries condition and diagnostic check up at start-up Settings password	•	lan touch lan touch	Ian, net, seal VP, poly VP	data data data data data	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	
Microprocessor Multilanguage menu Starting photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting Packaging programs setting Expiry dates setting Stand by Energy saving Safety thermostat On board battery ensures data retention even when the machine is switched off Peripheries condition and diagnostic check up at start-up Settings password Maintenance password	•	lan touch lan touch	Ian, net, seal VP, poly VP	data data data data data	Medical H D.A., SI Med/Vmed	
Microprocessor Multilanguage menu Starting photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting Packaging programs setting Expiry dates setting Stand by Energy saving Safety thermostat On board battery ensures data retention even when the machine is switched off Peripheries condition and diagnostic check up at start-up Settings password Maintenance password Display shutdown	•	lan touch lan touch	Ian, net, seal VP, poly VP	data data data data data seal, rebi	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	
Microprocessor Multilanguage menu Starting photocell Printing photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting Packaging programs setting Expiry dates setting Stand by Energy saving Safety thermostat On board battery ensures data retention even when the machine is switched off Peripheries condition and diagnostic check up at start-up Settings password Maintenance password Display shutdown Absolute counter	•	lan touch lan touch	Ian, net, seal VP, poly VP	data data data data data seal, rebi seal, rebi	Medical H D.A., SI Med/Vmed	
Microprocessor Multilanguage menu Starting photocell Printing photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting Packaging programs setting Expiry dates setting Stand by Energy saving Safety thermostat On board battery ensures data retention even when the machine is switched off Peripheries condition and diagnostic check up at start-up Settings password Maintenance password Display shutdown Absolute counter Settable counter	•	lan touch lan touch lan touch	Ian, net lan, net lan	data data data data data seal, rebi	Medical H D.A., SI Med/Vmed	
Microprocessor Multilanguage menu Starting photocell Printing photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting Packaging programs setting Expiry dates setting Stand by Energy saving Safety thermostat On board battery ensures data retention even when the machine is switched off Peripheries condition and diagnostic check up at start-up Settings password Maintenance password Display shutdown Absolute counter Settable counter Maintenance counter	•	lan touch lan touch lan touch	Ian, net, seal VP, poly VP	data data data data data seal, rebi seal, rebi seal	Medical H D.A., SI Med/Vmed	
Microprocessor Multilanguage menu Starting photocell Printing photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting Packaging programs setting Expiry dates setting Stand by Energy saving Safety thermostat On board battery ensures data retention even when the machine is switched off Peripheries condition and diagnostic check up at start-up Settings password Maintenance password Display shutdown Absolute counter Settable counter Maintenance counter Clock function		lan touch lan touch	Ian, net, seal VP, poly VP	data data data data data data seal, rebi seal, rebi seal data, rebi Prt	Medical H D.A., SI Med/Vmed Medical H D.A., SI Med/Vmed	
Microprocessor Multilanguage menu Starting photocell Printing photocell Printing photocell Speed adjustment Temperature setting Speed setting Time setting Operators setting Batchs setting Catalogues setting Packaging programs setting Expiry dates setting Stand by Energy saving Safety thermostat On board battery ensures data retention even when the machine is switched off Peripheries condition and diagnostic check up at start-up Settings password Maintenance password Display shutdown Absolute counter Settable counter Maintenance counter	•	lan touch lan touch lan touch	Ian, net lan, net lan	data data data data data seal, rebi seal, rebi seal	Medical H D.A., SI Med/Vmed	



SI Vmed

MODELS minirò H-lan minirò H-data Evo Medical H- 460/610 **DN 350 Vseal.it Pura** minirò H-lan touch minirò H-net Evo minirò H-seal Evo **Medical H 460/610 GD 301** Evo minirò H-seal VP Evo rebi Prt Evo Digital Accutemp minirò-poly VP rebi Evo SI Med / SI Vacuummed 520/720/1020 rebi VP Evo **VACUUM FUNCTION** Venturi vacuum pump SI Vmed Air suction capacity m3/h 11,3 SI Vmed Gas injection SI Vmed Vacuum time function SI Vmed Gas time function SI Vmed 9 settable programs

INTERFACES

INTENTAGES						
Display	4,3" colors	4,3" colors Pura	lan, net - 4 lines 20 cht	data, seal - 2 lines 16 cht	Medical H D.A., SI Med/Vmed	
		7" colors lan touch	seal VP, poly VP -2 lines 16 cht	rebi Prt, rebi VP - 2 lines 8 cht	4 lines 20 cht	
			rebi VP -2 lines 8 cht			
Touch screen	•	•				
Membrane keyboard			•	•	Medical H D.A., SI Med/Vmed	
Windows® CE operating system		lan touch				
PC connection	•	•	lan, net		Medical H D.A., SI Med/Vmed	
Remote control via VNC/Ethernet		lan touch				
Remote control					Medical H D.A., SI Med/Vmed	
Labels printer		•	lan, net, seal VP, poly VP			
Barcode scanner	•	Pura	net			
Lan net		•	lan			
RS 232 port	•	•	lan, net, seal VP, poly VP		Medical H D.A.	
USB port	•	•	net, seal VP, poly VP, rebi VP			
DIAG port					SI Med/Vmed	
Ethernet port		•	lan			
Configuration data for Ethernet connection		•	lan			
Back up recordings on PC		•			Medical H D.A., SI Med/Vmed	

BUILT-IN PRINTER FUNCTIONS

DUILI-IN PRINTER FUNCTIONS				
Built-in printer	•	lan, net, rebi VP,	data, rebi Prt	
On/off print	•	lan, net	data, rebi Prt	
EN symbols	•	lan, net	data, rebi Prt	
CE symbols	•	lan, net		
Print lines	Pura=1 line - lan touch=2 lines	lan=2 lines - net=1 line	data=1 line - rebi Prt =1 line	
Barcode	lan touch	lan		
Adjustable print position	•	lan, net	data, rebi Prt	
Standard or 180° rotation print image	•	lan, net		
Print automatic adjustment to pouch width	•	lan		
Settable printing formats	•	lan, net		
Settable character dimensions	•	lan, net	data	
Settable printing sequence	•	lan, net		
Production and expiry date	•	lan, net	data, rebi Prt	
Batch and operator alphanumeric code	•	lan, net	data	
Catalogue strings	•	lan, net		
Pouch progressive number	•	lan, net		
Pouch and sealer ID	•	lan, net		
Sealing speed	•	lan, net		
Sealing temperature	•	lan, net		
Sealing force	•	lan, net		
Sealing report print	•	lan, net		



MODELS	Vseal.it	Pura minirò H-lan touch	minirò H-lan minirò H-net Evo minirò H-seal VP Evo minirò-poly VP rebi VP Evo	minirò H-data Evo minirò H-seal Evo rebi Prt Evo rebi Evo	Medical H- 460/610 Medical H 460/610 Digital Accutemp SI Med / SI Vacuummed 520/720/1020	DN 350 GD 301 <i>Evo</i>
LABELS PRINTER FUNCTIONS		•	les set end VD selv VD			
Labels printer		•	lan, net, seal VP, poly VP			
On/off print		•	lan, net, seal VP, poly VP			
CE symbols		•	lan, net, seal VP, poly VP			
Barcode		•	lan, net, seal VP, poly VP			
Label dimensions setting			lan, net, seal VP, poly VP			
Standard or 180° rotation print image		lan touch	lan, seal VP, poly VP			
Multi labels print		lan touch	seal VP, poly VP lan, net, seal VP, poly VP			
Settable printing sequence		•				
Production and expiry date		•	lan, net, seal VP, poly VP lan, net, seal VP, poly VP			
Batch and operator alphanumeric code		•				
Catalogue strings		•	lan, net, seal VP, poly VP			
Pouch progressive number Pouch and sealer ID		•	lan, net, seal VP, poly VP lan, net, seal VP, poly VP			
		•				
Sealing speed		•	lan, net, seal VP, poly VP			
Sealing temperature		•	lan, net, seal VP, poly VP			
Sealing force		•	lan, net, seal VP, poly VP			
Sealing report print		•	lan, net, seal VP, poly VP			
Labels printer Barcode Scanner	•	•	lan, net, seal VP, poly VP net			
Traceability through Ethernet Traceability through serial port Traceability through local memory Traceability through USB pen drive	•	lan touch lan touch Pura	net, rebi VP, seal VP, poly VP		Medical H D.A. SI Med/Vmed	
Traceability through DIAG interface SECURSCAN Data Upload Software	•	Pura	4		Si ivieu/villeu	
VNC Remote Connection Software			net			
ROPEX Visualization Remote Connection Software		lan touch			Medical H D.A., SI Med/Vmed	
DATABASE						
Texts saving		•	lan, net, seal VP, rebi VP, poly VP			
Sealer ID saving		•	lan, net, seal VP, rebi VP, poly VP			
			ian, non sour vi, robi vi, pory vi			
<u> </u>		•	lan net seal VP rehi VP noly VP	data	The state of the s	
Operators saving		•	lan, net, seal VP, rebi VP, poly VP	data data_rehi Prt		
Operators saving Expiry date saving			lan, net, seal VP, poly VP	data, rebi Prt		
Operators saving Expiry date saving Batch code saving	•	•				
Operators saving Expiry date saving Batch code saving Presettable programs	•	•	lan, net, seal VP, poly VP	data, rebi Prt		
Operators saving Expiry date saving Batch code saving Presettable programs POWER SUPPLY		•	lan, net, seal VP, poly VP lan, net, seal VP, rebi VP,	data, rebi Prt data	100. 110/115. 220/240	DN 260 - 110/115 20
Operators saving Expiry date saving Batch code saving Presettable programs POWER SUPPLY Power supply (VAC) single phase	100; 110/115; 230/240	100; 110/115; 230/240	lan, net, seal VP, poly VP lan, net, seal VP, rebi VP, 100; 110/115; 230/240	data, rebi Prt data 100; 110/115; 230/240	100; 110/115; 230/240	DN 350 = 110/115-230 GD 301 <i>Evo</i> = 230
Operators saving Expiry date saving Batch code saving Presettable programs POWER SUPPLY Power supply (VAC) single phase Frequency (Hz)	100; 110/115; 230/240 50/60	100; 110/115; 230/240	lan, net, seal VP, poly VP lan, net, seal VP, rebi VP, 100; 110/115; 230/240 50/60	data, rebi Prt data 100; 110/115; 230/240 50/60	50/60	GD 301 <i>Evo</i> = 230 50/60
Operators saving Expiry date saving Batch code saving Presettable programs POWER SUPPLY Power supply (VAC) single phase Frequency (Hz) Multi-frequency	100; 110/115; 230/240 50/60	100; 110/115; 230/240 50/60	lan, net, seal VP, poly VP lan, net, seal VP, rebi VP, 100; 110/115; 230/240 50/60	data, rebi Prt data 100; 110/115; 230/240 50/60	50/60 •	GD 301 <i>Evo</i> = 230 50/60
Operators saving Expiry date saving Batch code saving Presettable programs POWER SUPPLY Power supply (VAC) single phase Frequency (Hz) Multi-frequency	100; 110/115; 230/240 50/60	100; 110/115; 230/240 50/60	lan, net, seal VP, poly VP lan, net, seal VP, rebi VP, 100; 110/115; 230/240 50/60 an = 650	data, rebi Prt data 100; 110/115; 230/240 50/60 data = 600	50/60	GD 301 <i>Evo</i> = 230 50/60 ● DN 350 = 900
Operators saving Expiry date saving Batch code saving Presettable programs POWER SUPPLY Power supply (VAC) single phase Frequency (Hz) Multi-frequency	100; 110/115; 230/240 50/60	100; 110/115; 230/240 50/60	lan, net, seal VP, poly VP lan, net, seal VP, rebi VP, 100; 110/115; 230/240 50/60 lan = 650 net = 600	data, rebi Prt data 100; 110/115; 230/240 50/60 • data = 600 seal = 600	50/60 ◆ Medical H = 1700/2100 Medical H D.A., = 2100/2500	GD 301 <i>Evo</i> = 230 50/60
Operators saving Expiry date saving Batch code saving Presettable programs	100; 110/115; 230/240 50/60	100; 110/115; 230/240 50/60	lan, net, seal VP, poly VP lan, net, seal VP, rebi VP, 100; 110/115; 230/240 50/60 an = 650	data, rebi Prt data 100; 110/115; 230/240 50/60 data = 600	50/60	GD 301 <i>Evo</i> = 230 50/60 • DN 350 = 900



MODELS	Vseal.it	Pura minirò H-lan touch	minirò H-lan minirò H-net Evo minirò H-seal VP Evo minirò-poly VP rebi VP Evo	minirò H-data Evo minirò H-seal Evo rebi Prt Evo rebi Evo	Medical H- 460/610 Medical H 460/610 Digital Accutemp SI Med / SI Vacuummed 520/720/1020	DN 350 GD 301 <i>Evo</i>
ACCESSORIES	•	•	•	•	•	DNIOCO
Front Support Tray	•	•	•	•	•	DN 350
Roller Support Tray	•	•	•	•	•	•
Roll holder		•	_	•	•	•
Motorized transport belt		•	poly VP			
Thermal transfer labels printer			lan, net, seal VP, poly VP			
Barcode Scanner	•	•	net			
Ethernet module		Pura	lan			
Securscan Data Upload Software	•	Pura	net			
SERVICES	•	•	•	•	•	
Seal strength test EN 868-5			-	-	-	•
Machine calibration	•	•	•	•	Medical H D.A., SI Med/Vmed	
			-	-	Medical H D.A., SI Med/Vmed	•
Machine calibration	•	•	•	•	Medical H D.A., SI Med/Vmed	
Machine calibration Maintenance DIMENSIONS / WEIGHTS	•	•	•	•	Medical H D.A., SI Med/Vmed	
Machine calibration Maintenance DIMENSIONS / WEIGHTS	•	•	•	•	Medical H D.A., SI Med/Vmed	•
Machine calibration Maintenance DIMENSIONS / WEIGHTS	•	• • Pura 553x264x192	lan 621x300x330	data 554x270x222	Medical H D.A., SI Med/Vmed Medical H 460 /Medical H	DN 439x311x175
Machine calibration Maintenance DIMENSIONS / WEIGHTS	•	• • Pura 553x264x192	lan 621x300x330 net 554x270x222	data 554x270x222 seal 554x270x222	Medical H D.A., SI Med/Vmed Medical H 460 /Medical H D.A 460 560x360x220	DN 439x311x175
Machine calibration Maintenance DIMENSIONS / WEIGHTS	•	• • Pura 553x264x192	lan 621x300x330 net 554x270x222 seal VP 554x270x222	data 554x270x222 seal 554x270x222 rebi Prt 473x235x181	Medical H D.A., SI Med/Vmed Medical H 460 /Medical H D.A 460 560x360x220 Medical H 610 /Medical H	DN 439x311x175
Machine calibration Maintenance DIMENSIONS / WEIGHTS	•	• • Pura 553x264x192	lan 621x300x330 net 554x270x222 seal VP 554x270x222 rebi VP 473x235x181	data 554x270x222 seal 554x270x222 rebi Prt 473x235x181	Medical H D.A., SI Med/Vmed Medical H 460 /Medical H D.A 460 560x360x220 Medical H 610 /Medical H D.A 610 710x360x220	DN 439x311x175
Machine calibration Maintenance DIMENSIONS / WEIGHTS Dimensions (LxDxH mm)	•	• • Pura 553x264x192	lan 621x300x330 net 554x270x222 seal VP 554x270x222 rebi VP 473x235x181	data 554x270x222 seal 554x270x222 rebi Prt 473x235x181	Medical H D.A., SI Med/Vmed Medical H 460 /Medical H D.A 460 560x360x220 Medical H 610 /Medical H D.A 610 710x360x220 SI Med/Vmed 590x487x241 -	DN 439x311x175
Machine calibration Maintenance DIMENSIONS / WEIGHTS Dimensions (LxDxH mm)	590x520x300	Pura 553x264x192 lan touch 621x300x330	lan 621x300x330 net 554x270x222 seal VP 554x270x222 rebi VP 473x235x181 poly VP 554x270x237	data 554x270x222 seal 554x270x222 rebi Prt 473x235x181 rebi 473x235x181	Medical H D.A., SI Med/Vmed Medical H 460 /Medical H D.A 460 560x360x220 Medical H 610 /Medical H D.A 610 710x360x220 SI Med/Vmed 590x487x241 790x487x241 - 1110x487x232	DN 439x311x175 GD 301 <i>Evo</i> 452x412x307
Machine calibration Maintenance DIMENSIONS / WEIGHTS Dimensions (LxDxH mm)	590x520x300	Pura 553x264x192 lan touch 621x300x330 Pura 18 Kg	lan 621x300x330 net 554x270x222 seal VP 554x270x222 rebi VP 473x235x181 poly VP 554x270x237	data 554x270x222 seal 554x270x222 rebi Prt 473x235x181 rebi 473x235x181	Medical H D.A., SI Med/Vmed Medical H 460 /Medical H D.A 460 560x360x220 Medical H 610 /Medical H D.A 610 710x360x220 SI Med/Vmed 590x487x241 - 790x487x241 - 1110x487x232 Medical H 30/37 Kg	DN 439x311x175 GD 301 <i>Evo</i> 452x412x307
Machine calibration Maintenance	590x520x300	Pura 553x264x192 lan touch 621x300x330 Pura 18 Kg	lan 621x300x330 net 554x270x222 seal VP 554x270x222 rebi VP 473x235x181 poly VP 554x270x237 lan 26 Kg net 18,5 Kg	data 554x270x222 seal 554x270x222 rebi Prt 473x235x181 rebi 473x235x181 data 18 Kg seal 15 Kg	Medical H D.A., SI Med/Vmed Medical H 460 /Medical H D.A 460 560x360x220 Medical H 610 /Medical H D.A 610 710x360x220 SI Med/Vmed 590x487x241 - 790x487x241 - 1110x487x232 Medical H 30/37 Kg Medical H D.A. 30/37 Kg	DN 439x311x175 GD 301 <i>Evo</i> 452x412x307



For more then 60 years Gandus Saldatrici has been a point of reference for the production of heat sealers and packaging lines

for packaging in heat sealable bags and sacks. Its steady development is proven by its successful presence in the most technologically advanced packaging markets: medical, vacuum, modified atmosphere, dosing systems, automatic packaging machines, complete lines.

Today Gandus Saldatrici can offer you a wide and complete range of products, which includes 33 different models available in more than 100 versions.

With its internationally developed sales network, Gandus Saldatrici is now represented in more than 30 countries all over the world

Machines are designed and built under the most strict qualitative controls and in accordance with ISO 9001:2015 regulations.

Quality is our vocation.







Medical Sealers Division

Gandus Saldatrici S.r.l. socio unico Via Milano, 5 - 20010 Cornaredo (Milano) - Italy Tel. +39 02 93.19.41 - Fax + 39 02 93.56.88.03 info@gandus.it - www.gandus.it



